

DESCRIPTION


Species Reactivity	Human/Mouse
Specificity	Detects human and mouse Wnt-8b in direct ELISAs and Western blots. In direct ELISAs, approximately 25% cross-reactivity with recombinant mouse Wnt-8a is observed.
Source	Polyclonal Goat IgG
Purification	Antigen Affinity-purified
Immunogen	<i>E. coli</i> -derived recombinant mouse Wnt-8b Trp22-Glu52, Leu195-Arg269 Accession # Q9WUD6
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. See Certificate of Analysis for details. *Small pack size (-SP) is supplied as a 0.2 µm filtered solution in PBS.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the [Technical Information](#) section on our website.

	Recommended Concentration	Sample
Western Blot	1 µg/mL	See Below

DATA



PREPARATION AND STORAGE

Reconstitution	Reconstitute at 0.2 mg/mL in sterile PBS.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below. *Small pack size (-SP) is shipped with polar packs. Upon receipt, store it immediately at -20 to -70 °C
Stability & Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> ● 12 months from date of receipt, -20 to -70 °C as supplied. ● 1 month, 2 to 8 °C under sterile conditions after reconstitution. ● 6 months, -20 to -70 °C under sterile conditions after reconstitution.

BACKGROUND

Mouse Wnt-8b is a presumably 40-50 kDa secreted glycoprotein member of the Wnt family of proteins. Mature mouse Wnt-8b is 329 amino acids (aa) in length and contains 24 cysteines plus three potential N-linked glycosylation sites. Little is known about Wnt-8b. It appears to play a role in forebrain development. The mature mouse protein is 100%, 98% and 97% aa identical to mature rat, canine and human Wnt-8b, respectively.