Rat IL-18/IL-1F4 Antibody

Antigen Affinity-purified Polyclonal Goat IgG Catalog Number: AF521


| DESCRIPTION | | | |
|--------------------|---|--|--|
| Species Reactivity | Rat | | |
| Specificity | Detects rat IL-18/IL-1F4 in direct ELISAs and Western blots. In Western blots, approximately 35% cross-reactivity with recombinant mouse IL-18 is observed and approximately 10% cross-reactivity with recombinant human IL-18 is observed. | | |
| Source | Polyclonal Goat IgG | | |
| Purification | Antigen Affinity-purified | | |
| Immunogen | E. coli-derived recombinant rat IL-18/IL-1F4 His37-Ser194 Accession # P97636 | | |
| Endotoxin Level | <0.10 EU per 1 µg of the antibody by the LAL method. | | |
| Formulation | Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. See Certificate of Analysis for details. *Small pack size (-SP) is supplied as a 0.2 µm filtered solution in PBS. | | |

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. General Protocols are available in the Technical Information section on our website.

| | Recommended Concentration | Sample |
|----------------|------------------------------|---|
| Western Blot | 0.1 μg/mL | Recombinant Rat IL-18/IL-1F4 (Catalog # 521-RL) |
| Neutralization | • | y to neutralize IL-18/IL-1F4-induced IFN- γ secretion in activated mouse T cells [Ahn, H.J. et . 159 :2125]. The Neutralization Dose (ND ₅₀) is typically 1-3 μ g/mL in the presence of |
| | 15 ng/mL Recombina | nt Rat IL-18/IL-1F4 and 0.1 ng/mL Recombinant Mouse IL-12. |

DATA



IFN-γ Secretion Induced by IL-18/IL-1F4 and Neutralization by Rat IL-18/IL-1F4 Antibody. In the presence of Recombinant Mouse IL-12 (0.1 ng/mL, Catalog # 419-ML), Recombinant Rat IL-18/IL-1F4 (Catalog # 521-RL) stimulates IFN-γ secretion in activated mouse T cells in a dose-dependent manner (orange line), as measured by the Mouse IFN-γ Quantikine ELISA Kit (Catalog # MIF00). Under these conditions, IFN-γ secretion elicited by Recombinant Rat IL-18/IL-1 F 4 (15 ng/mL) is neutralized (green line) by increasing concentrations of Rat IL-18/IL-1F4 Antigen Affinitypurified Polyclonal Antibody (Catalog # AF521). The ND_{50} is typically 1-3 µg/mL.

| PREPARATION AND STORAGE | | |
|-------------------------|--|--|
| Reconstitution | Reconstitute at 0.2 mg/mL in sterile PBS. | |
| Shipping | The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below. *Small pack size (-SP) is shipped with polar packs. Upon receipt, store it immediately at -20 to -70 °C | |
| Stability & Storage | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. 12 months from date of receipt, -20 to -70 °C as supplied. 1 month, 2 to 8 °C under sterile conditions after reconstitution. 6 months, -20 to -70 °C under sterile conditions after reconstitution. | |





Rat IL-18/IL-1F4 Antibody

Antigen Affinity-purified Polyclonal Goat IgG Catalog Number: AF521

BACKGROUND

Interleukin-18 (IL-18), also known as IL-1F4 and IFN-γ inducing factor (IGIF), is a member of the IL-1 family of cytokines and is a key molecule in the innate immune response (1). Rat IL-18 is synthesized as a 24 kDa proprotein that contains a 36 amino acid (aa) propeptide and a 158 aa mature region (2). Under inflammatory conditions, the propeptide is cleaved by Caspase-1 in the cytoplasm to liberate the mature nonglycosylated 18 kDa monomeric IL-18 (3, 4). Mature rat IL-18 shares 91% aa sequence identity with mouse IL-18 and 60-64% aa sequence identity with human, canine, feline, porcine, and rhesus macaque IL-18. IL-18 is secreted by a variety of cell types including macrophages, dendritic cells, and epithelial cells (1, 5). Circulating mature IL-18 is sequestered by soluble IL-18 binding proteins (IL-18 BP) that inhibit IL-18 bioactivity (6). IL-18 interacts with the widely expressed IL-18 Rα which then recruits the signaling subunit IL-18 Rβ (7, 8). The IL-1 family member IL-1F7 also binds to IL-18 Rα but does not recruit IL-18 Rβ or induce signaling (9). IL-1F7 binds IL-18 BP and enhances its neutralizing effect on IL-18 activity (9). IL-18 synergizes with other cytokines to activate NK, Th1, and Th17 cells and to increase the production of IFN-γ (1, 5, 10-12). IL-18 can also promote Th2 cytokine release which reduces the effectiveness of antiviral responses (13, 14). Increased levels of active IL-18 stimulates Th1 and NK cells to target tumor cells, but it can also promote angiogenesis, metastasis, and tumor cell immune evasion (11).

References:

- 1. Arend, W.P. et al. (2008) Immunol. Rev. 223:20.
- 2. Culhane, A.C. et al. (1998) Mol. Psychiatry 3:362.
- 3. Ghayur, T. et al. (1997) Nature 386:619.
- 4. Gu, Y. et al. (1997) Science 275:206.
- 5. Boraschi, D. and C.A. Dinarello (2006) Eur. Cytokine Netw. 17:224.
- 6. Novick, D. et al. (1999) Immunity 10:127.
- 7. Torigoe, K. et al. (1997) J. Biol. Chem. 272:25737.
- 8. Born, T.L. et al. (1998) J. Biol. Chem. 273:29445.
- 9. Bufler, P. et al. (2002) Proc. Natl. Acad. Sci. USA 99:13723.
- 10. Takeda, K. et al. (1998) Immunity 8:383.
- 11. Park, S. et al. (2007) Cell. Mol. Immunol. 4:329.
- 12. Yoshimoto, T. et al. (1998) J. Immunol. 161:3400.
- 13. Hoshino, T. et al. (2001) J. Immunol. 166:7014.
- 14. Iannello, A. et al. (2009) AIDS Rev. 11:115.
- 15. Rabkin, S.W. (2009) Nat. Clin. Pract. Cardiovasc. Med. 6:192.
- 16. Netea, M.G. et al. (2006) Nat. Med. 12:650.

