

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human Mindin in Western blots.
Source	Polyclonal Goat IgG
Purification	Antigen Affinity-purified
Immunogen	Mouse myeloma cell line NS0-derived recombinant human Mindin Glu32-Val331 Accession # Q9BUD6
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS with BSA as a carrier protein. See Certificate of Analysis for details.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

	Recommended Concentration	Sample
Western Blot	0.1 µg/mL	Recombinant Human Mindin

PREPARATION AND STORAGE

Reconstitution	Reconstitute at 0.2 mg/mL in sterile PBS.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	<p>Use a manual defrost freezer and avoid repeated freeze-thaw cycles.</p> <ul style="list-style-type: none"> ● 12 months from date of receipt, -20 to -70 °C as supplied. ● 1 month, 2 to 8 °C under sterile conditions after reconstitution. ● 6 months, -20 to -70 °C under sterile conditions after reconstitution.

BACKGROUND

MindIn, also known as Spondin-2 and DIL-1 (Differentially expressed in cancerous and non-cancerous lung cells 1) belongs to the F-spondin family of secreted extracellular matrix proteins. Family members are characterized by two N-terminal F-spondin domains (FS1 and FS2), and a third C-terminal thrombospondin type 1 repeat (TSP-1). Secreted Mindin exists as a disulfide-linked dimer that also forms higher order oligomers. It functions as a pattern recognition receptor for microbes and an adhesion molecule for neurons. Mindin is secreted by neurons, macrophages and fibroblasts in a variety of tissues. Human Mindin shares 88% amino acid sequence homology with mouse and rat proteins.