



***Affinity Purified
Donkey Anti-sheep IgG (H & L) Antibody***

Catalog Number: D-501-C-ABS

Lot Numbers: WBM0308091

WBM0312091

Product Specifications

Amount	◆ 0.5 mg
Form	◆ Lyophilized
Purity	◆ Affinity purified
Buffer	◆ 10 mM sodium phosphate, 0.85% (w/v) sodium chloride, pH 7.2 (PBS) with 5% trehalose
Preservative	◆ 0.05% (w/v) sodium azide
Reconstitution	◆ Reconstitute with sterile PBS. If 0.5 mL of PBS is used, the antibody concentration will be 1 mg/mL.
Storage	◆ 2° - 8° C
Specificity	◆ This antibody reacts with heavy chains on sheep IgG, and with light chains on all sheep immunoglobulins. No reactivity is observed to non-immunoglobulin sheep serum proteins. This antibody is adsorbed against solid phase human, mouse, bovine, and rabbit serum and human and rabbit IgG.
Note	◆ Sodium azide may react with lead and copper plumbing to form explosive metal azides. Flush with large amounts of water during disposal.