

DESCRIPTION


Species Reactivity	Human
Specificity	Recognizes human Leukotriene B4 R1 (Pettersson, A. <i>et al.</i> (2000) <i>Biochem. Biophys. Res. Commun.</i> 279 :520).
Source	Monoclonal Mouse IgG ₁ Clone # 203/14F11
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	HeLa cervical epithelial carcinoma cell line transfected with human Leukotriene B4 R1
Conjugate	Alexa Fluor 700 Excitation Wavelength: 675-700 nm Emission Wavelength: 723 nm
Formulation	Supplied at 0.2 mg/mL in a saline solution containing BSA and Sodium Azide. See Certificate of Analysis for details. *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

	Recommended Concentration	Sample
Flow Cytometry	0.25-1 µg/10 ⁶ cells	See Below

DATA



PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. ● 12 months from date of receipt, 2 to 8 °C as supplied.

BACKGROUND

Polymorphonuclear granulocytes secrete the lipid chemotactic mediator Leukotriene B4 (LTB₄) in response to inflammatory stimuli (1). Neutrophils, monocytes and lymphocytes respond to LTB₄ via specific receptors localized on the cell surface (2-4). The high affinity LTB₄ Receptor known as BLT1 is only expressed on leukocytes (5-7) while a second low affinity receptor BLT2 is expressed more ubiquitously (8, 9). The BLT1 and BLT2 are G-protein linked seven-transmembrane spanning receptors that share about 37-45% amino acid identity (8, 9).

References:

- Samuelsson, B. *et al.* (1987) *Science* **237**:1171.
- Brom, J. and W. Konig (1989) *Immunology* **68**:479.
- Patry, C. *et al.* (1996) *Prostaglandins, Leukotrienes and Essential Fatty Acids* **54**:361.
- Dasari, V.R. *et al.* (2000) *Immunopharm.* **48**:157.
- Pettersson, A. *et al.* (2000) *Biochem. Biophys. Res. Comm.* **279**:520.
- Yokomizo, T. *et al.* (1997) *Nature* **387**:620.
- Kato, K. *et al.* (2000) *J. Exp. Med.* **192**:413.
- Yokomizo, T. *et al.* (2000) *J. Exp. Med.* **192**:421.
- Kamohara, M. *et al.* (2000) *J. Biol. Chem.* **275**:27000.

Human Leukotriene B4 R1 Alexa Fluor® 700-conjugated Antibody

Monoclonal Mouse IgG₁ Clone # 203/14F11

Catalog Number: FAB099N

100 µg

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc. and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.