

## DESCRIPTION


<b>Species Reactivity</b>	Human/Mouse
<b>Specificity</b>	Recognizes a carbohydrate epitope of SSEA-4 (1, 2).
<b>Source</b>	Monoclonal Mouse IgG <sub>3</sub> Clone # MC-813-70
<b>Purification</b>	Protein A or G purified from hybridoma culture supernatant
<b>Immunogen</b>	2120Ep human embryonal carcinoma cell line
<b>Conjugate</b>	Allophycocyanin Excitation Wavelength: 620-650 nm Emission Wavelength: 660-670 nm
<b>Formulation</b>	Supplied in a saline solution containing BSA and Sodium Azide. See Certificate of Analysis for details.  *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

## APPLICATIONS

**Please Note:** Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

	Recommended Concentration	Sample
<b>Flow Cytometry</b>	10 µL/10 <sup>6</sup> cells	See Below

## DATA



## PREPARATION AND STORAGE

**Shipping** The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.

**Stability & Storage** **Protect from light. Do not freeze.**

- 12 months from date of receipt, 2 to 8 °C as supplied.

## BACKGROUND

SSEA-4 is expressed on the surface of human embryonal carcinoma (EC) cells (the pluripotent stem cells of teratocarcinomas), human embryonic germ cells (EG), and human embryonic stem cells (ES). Expression of SSEA-4 is down-regulated following differentiation of human EC cells. In contrast, the differentiation of murine EC and ES cells may be accompanied by an increase in SSEA-4 expression (1-4).

### References:

1. Shevinsky, L.H. *et al.* (1982) *Cell* **30**:697.
2. Kannagi, R. *et al.* (1983) *EMBO J.* **2**:2355.
3. Thomson, J.A. and J.S. Odorico (2000) *Trends Biotechnol.* **18**:53.
4. Draper, J.S. *et al.* (2002) *J. Anat.* **200**:249.