

DESCRIPTION


Species Reactivity	Human
Specificity	Detects human TRA-1-85 antigen.
Source	Monoclonal Mouse IgG ₁ Clone # TRA-1-85
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	2120Ep human embryonal carcinoma cell line
Conjugate	Phycoerythrin Excitation Wavelength: 488 nm Emission Wavelength: 565-605 nm
Formulation	Supplied in a saline solution containing BSA and Sodium Azide. See Certificate of Analysis for details. *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

	Recommended Concentration	Sample
Flow Cytometry	10 µL/10 ⁶ cells	See Below

DATA



PREPARATION AND STORAGE

Shipping The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.

Stability & Storage **Protect from light. Do not freeze.**

- 12 months from date of receipt, 2 to 8 °C as supplied.

BACKGROUND

The TRA-1-85 antigen, also known as OKa blood group antigen, is a specific epitope within the protein known as Basigin, EMMPRIN and CD147. It is a cell surface determinant expressed on almost all human cell types. This antibody has been used in somatic cell hybrid studies to identify tissues of partial human origin (1, 2).

References:

1. Williams, B.P. *et al.* (1988) *Immunogenetics*. **27**:322.
2. Spring, F.A. *et al.* (1997) *Eur. J. Immunol.* **27**:891.