

Human S1P₅/EDG-8 Alexa Fluor® 647-conjugated Antibody

Monoclonal Mouse IgG_{2B} Clone # 282503

Catalog Number: FAB3964R

100 µg

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human S1P ₅ /EDG-8. Stains human S1P ₅ /EDG-8 transfectants but not irrelevant transfectants.
Source	Monoclonal Mouse IgG _{2B} Clone # 282503
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	BaF3 mouse pro-B cell line transfected with human S1P ₅ /EDG-8 Met1-Asp398 Accession # Q9H228
Conjugate	Alexa Fluor 647 Excitation Wavelength: 650 nm Emission Wavelength: 668 nm
Formulation	Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide. See Certificate of Analysis for details. *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

	Recommended Concentration	Sample
Flow Cytometry	0.25-1 µg/10 ⁶ cells	HEK293 human embryonic kidney cell line transfected with human S1P ₅ /EDG-8 and eGFP

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Protect from light. Do not freeze. <ul style="list-style-type: none"> 12 months from date of receipt, 2 to 8 °C as supplied.

BACKGROUND

S1P₅ is also known as EDG-8 and nerve growth factor-related G-protein-coupled receptor-1 (NRG-1). S1P₅ is a 398 amino acid (aa) seven-transmembrane receptor putative glycoprotein that binds the lysolipid phosphoric acid mediator, sphingosine 1-phosphate. Extracellular portions of human S1P₅ show 96% and 97% aa identity with mouse and rat S1P₅, respectively. Isoform 1 is expressed at a low level in peripheral tissues. Isoform 2 has an alternate C-terminal that is 88 aa shorter and is expressed mainly in brain, spleen, and PBMC. S1P₅ is upregulated in large granular lymphocytic leukemias.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.