

Monoclonal Anti-human CD1a-Phycoerythrin

Catalog Number: FAB7076P

Lot Number: ACQL02

100 Tests

Reagents Provided

Phycoerythrin (PE)-conjugated mouse monoclonal anti-human CD1a: Supplied as 10 µg of antibody in 1 mL saline containing up to 0.5% BSA and 0.1% sodium azide.

Clone #: 703217

Isotype: mouse IgG₁

Reagents Not Provided

- Flow Cytometry Staining Buffer (Catalog # FC001) or other BSA-supplemented saline buffer.

Storage


Reagents are stable for **twelve months** from the date of receipt when stored in the dark at 2-8 °C.

Intended Use

Designed to quantitatively determine the percentage of cells bearing CD1a within a population and qualitatively determine the density of CD1a on cell surfaces by flow cytometry.

Product Description

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified NS0-derived recombinant human CD1a (rhCD1a; aa 1-300; Accession # P06126) extracellular domain. The IgG fraction of the tissue culture supernatant was purified by Protein A or G affinity chromatography. The purified antibody was then conjugated to PE fluorochrome. Cell surface expression of CD1a is determined by flow cytometry using 488 nm wavelength excitation and monitoring emitted fluorescence with a detector optimized to collect peak emissions at 565-605 nm.



Molt-4 cells were stained with PE-conjugated anti-human CD1a (Catalog # FAB7076P, filled histogram) or isotype control (Catalog # IC002P, open histogram).

Background Information

CD1a is a 49 kDa transmembrane glycoprotein in the CD1 family of glycolipid antigen-presenting MHC-like molecules. CD1a contains one Ig-like domain in its extracellular region. It is expressed by most non-human mammals but not by mice or rats. Complexes of CD1a with β2-microglobulin and endogenous glycolipids are constitutively expressed on antigen-presenting cells, cortical thymocytes, and Langerhans cells. CD1a is a target of autoreactive Th22 helper T cells in the skin.

Flow Cytometry Validation

This antibody has been tested for flow cytometry using Molt-4 cells.

- Cells may be Fc-blocked with 1 µg of human IgG/10⁵ cells for 15 minutes at room temperature. Do not wash excess blocking IgG from this reaction.
- After blocking, 10 µL of conjugated antibody was added to up to 1 x 10⁶ cells and incubated for 30 minutes at room temperature.
- Unbound antibody was removed by washing the cells twice in Flow Cytometry Staining Buffer (Catalog # FC001). Note that whole blood requires a RBC lysis step at this point using Flow Cytometry Human Lyse Buffer (Catalog # FC002).
- The cells were resuspended in Flow Cytometry Staining Buffer for final flow cytometric analysis. As a control for this analysis, cells in a separate tube should be treated with PE-labeled mouse IgG₁ antibody. This procedure may need to be modified, depending upon the cell type and final utilization. Individual users may need to titrate to determine the optimal reagent amount for their specific use.

Warning: Contains sodium azide as a preservative - sodium azide may react with lead and copper plumbing to form explosive metal azides. Flush with large volumes of water during disposal.

FOR RESEARCH USE ONLY. NOT FOR USE IN HUMANS.

R&D Systems Inc.
1-800-343-7475