

DESCRIPTION

Species Reactivity	Rat
Specificity	Detects rat IgG ₁ , IgG _{2A} , IgG _{2B} heavy and light chains in direct ELISAs and Western blots. In direct ELISAs, less than 1% cross-reactivity with human IgG and rabbit IgG is observed.
Source	Polyclonal Goat IgG
Purification	Antigen Affinity-purified
Immunogen	Rat IgG
Formulation	Supplied as a 0.2 um filtered solution in Sodium Phosphate and Proclin 300. See Certificate of Analysis for details.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

Western Blot Use at a 1:1000 dilution.

PREPARATION AND STORAGE

Shipping The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.

Stability & Storage **Do not freeze.**

- 6 months from date of receipt, 2 to 8 °C as supplied.