

Sheep IgG Horseradish Peroxidase-conjugated Antibody

Antigen Affinity-purified Polyclonal Donkey IgG
Catalog Number: HAF016

DESCRIPTION

Species Reactivity	Sheep
Specificity	Detects sheep IgG heavy and light chains in direct ELISAs and Western blots. In Western blots, less than 5% cross-reactivity with mouse IgG, rabbit IgG, human IgG, and chicken IgY is observed.
Source	Polyclonal Donkey IgG
Purification	Antigen Affinity-purified
Immunogen	Sheep IgG
Formulation	Sodium Phosphate and Proclin 300. See Certificate of Analysis for details.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. [General Protocols](#) are available in the Technical Information section on our website.

Western Blot Use at a 1:1000 dilution.

PREPARATION AND STORAGE

Shipping	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	<p>Do not freeze.</p> <ul style="list-style-type: none"> 6 months from date of receipt, 2 to 8 °C as supplied.