

## DESCRIPTION


<b>Species Reactivity</b>	Human
<b>Specificity</b>	Detects human BLIMP1/PRDM1 in Western blots.
<b>Source</b>	Monoclonal Mouse IgG <sub>1</sub> Clone # 646702
<b>Purification</b>	Protein A or G purified from hybridoma culture supernatant
<b>Immunogen</b>	<i>E. coli</i> -derived recombinant human BLIMP1/PRDM1 Lys667-Cys789 Accession # O75626
<b>Conjugate</b>	Allophycocyanin Excitation Wavelength: 620-650 nm Emission Wavelength: 660-670 nm
<b>Formulation</b>	Supplied in a saline solution containing BSA and Sodium Azide. See Certificate of Analysis for details.  *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions.

## APPLICATIONS

**Please Note:** Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

	Recommended Concentration	Sample
Intracellular Staining by Flow Cytometry	10 μL/10 <sup>6</sup> cells	See Below

## DATA



## PREPARATION AND STORAGE

<b>Shipping</b>	The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below.
<b>Stability &amp; Storage</b>	<b>Protect from light. Do not freeze.</b> <ul style="list-style-type: none"> <li>● 12 months from date of receipt, 2 to 8 °C as supplied.</li> </ul>

## BACKGROUND

Human BLIMP1 (B lymphocyte-induced maturation protein 1, also known as PRDM1) is a 91 kDa zinc-finger transcriptional repressor that promotes B cell maturation into plasma cells. It is 789 amino acids (aa) in length, and contains an N-terminal S-E-T domain (aa 64-170) and four C-terminal C2H2-type zinc-finger motifs (aa 539-645). The SET domain interacts with chromatin modifiers, while the zinc fingers bind to DNA. There is one 80 kDa, 691 aa alternate splice form that utilizes an internal start site. This results in a substitution of three aa for the first 101 aa of the long form, and the loss of the SET domain. At least 10 mutations exist, resulting in proteins of 61-603 aa in length. Over aa 667-789, human BLIMP1 is 89% aa identical to mouse BLIMP1.