

Human/Mouse RUNX3/CBFA3 Alexa Fluor® 405-conjugated Antibody

Monoclonal Mouse IgG_{2A} Clone # 527327

Catalog Number: IC3765V
100 µg

DESCRIPTION

| | |
|---------------------------|--|
| Species Reactivity | Human/Mouse |
| Specificity | Detects human and mouse RUNX3/CBFA3 in Western blots. |
| Source | Monoclonal Mouse IgG _{2A} Clone # 527327 |
| Purification | Protein A or G purified from hybridoma culture supernatant |
| Immunogen | <i>E. coli</i> -derived recombinant human RUNX3/CBFA3 Lys186-Tyr415 Accession # Q13761 |
| Conjugate | Alexa Fluor 405 Excitation Wavelength: 405 nm Emission Wavelength: 421 nm |
| Formulation | Supplied 0.2 mg/mL in a saline solution containing BSA and Sodium Azide. See Certificate of Analysis for details. *Contains <0.1% Sodium Azide, which is not hazardous at this concentration according to GHS classifications. Refer to the Safety Data Sheet (SDS) for additional information and handling instructions. |

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

| | Recommended Concentration | Sample |
|---|---------------------------------|--|
| Intracellular Staining by Flow Cytometry | 0.25-1 µg/10 ⁶ cells | Human PBMC treated with PMA and mouse splenocytes fixed with paraformaldehyde and permeabilized with saponin |

PREPARATION AND STORAGE

| | |
|--------------------------------|--|
| Shipping | The product is shipped with polar packs. Upon receipt, store it immediately at the temperature recommended below. |
| Stability & Storage | Protect from light. Do not freeze. <ul style="list-style-type: none"> 12 months from date of receipt, 2 to 8 °C as supplied. |

BACKGROUND

RUNX3, also called CBFA3, AML-2 or PEBP2-αC, is a member of the Runt domain family of nuclear transcriptional regulators. All of the RUNX proteins form dimers with CBF-β. The runt domain (aa 54-186) is required for DNA binding, while a pro/ser/thr-rich region (aa 191-415) transcriptionally activates target genes. Isoform 2 has an alternate 19 aa in place of the N-terminal 5 aa of isoform 1. The 415 aa Human RUNX3 shares 91% aa identity with mouse or rat RUNX3. RUNX3 is necessary for growth control of gastric epithelium, neurogenesis of dorsal root ganglia, and T cell differentiation. RUNX3 expression is frequently mutated in tumors and appears to be silenced by methylation.

PRODUCT SPECIFIC NOTICES

This product is provided under an agreement between Life Technologies Corporation and R&D Systems, Inc, and the manufacture, use, sale or import of this product is subject to one or more US patents and corresponding non-US equivalents, owned by Life Technologies Corporation and its affiliates. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components (1) in manufacturing; (2) to provide a service, information, or data to an unaffiliated third party for payment; (3) for therapeutic, diagnostic or prophylactic purposes; (4) to resell, sell, or otherwise transfer this product or its components to any third party, or for any other commercial purpose. Life Technologies Corporation will not assert a claim against the buyer of the infringement of the above patents based on the manufacture, use or sale of a commercial product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.