

DESCRIPTION

| | |
|---------------------------|--|
| Species Reactivity | Mouse |
| Specificity | Detects mouse Chemerin in ELISAs. In sandwich immunoassays, no cross-reactivity or interference with recombinant human Chemerin, recombinant mouse (rm) Cystatin C, or rmFetuin A is observed. |
| Source | Monoclonal Rat IgG _{2B} Clone # 372409 |
| Purification | Protein A or G purified from hybridoma culture supernatant |
| Immunogen | <i>E. coli</i> -derived recombinant mouse Chemerin Thr17-Ser156 Accession # Q9DD06 |
| Formulation | Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. See Certificate of Analysis for details. *Small pack size (-SP) is supplied as a 0.2 µm filtered solution in PBS. |

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

| Mouse Chemerin Sandwich Immunoassay | | Reagent |
|-------------------------------------|---------------|--|
| ELISA Capture | 2-8 µg/mL | Mouse Chemerin Antibody (Catalog # MAB23251) |
| ELISA Detection | 0.5-2.0 µg/mL | Mouse Chemerin Biotinylated Antibody (Catalog # BAM2325) |
| Standard | | Recombinant Mouse Chemerin (Catalog # 2325-CM) |

PREPARATION AND STORAGE

| | |
|--------------------------------|--|
| Reconstitution | Reconstitute at 0.5 mg/mL in sterile PBS. |
| Shipping | The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below. *Small pack size (-SP) is shipped with polar packs. Upon receipt, store it immediately at -20 to -70 °C |
| Stability & Storage | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> ● 12 months from date of receipt, -20 to -70 °C as supplied. ● 1 month, 2 to 8 °C under sterile conditions after reconstitution. ● 6 months, -20 to -70 °C under sterile conditions after reconstitution. |

BACKGROUND

Mouse Chemerin, also known as Tazarotene-induced Gene-2 (TIG2), is a new, but distant member of the cystatin superfamily (1-3). Members of this superfamily contain at least two intrachain disulfide bonds and an α -helical structure over a distance of about 100 amino acids (aa) (2, 3). Chemerin is synthesized as a 162 aa precursor that contains a hydrophobic N-terminal sequence, an intervening 140 aa cystatin-fold containing domain, and a six aa C-terminal prosegment (4-6). Within the cystatin-fold domain there are three intrachain disulfide bonds that contribute to the characteristic fold (4, 7). The precursor molecule is described as undergoing proteolytic processing at both termini by unknown proteases. The N-terminal 16 residue hydrophobic segment is described as being either a signal sequence or a transmembrane (TM) segment for a type II TM protein (5, 8). In either case it gives rise to a soluble proform that undergoes further processing at the C-terminus (5). In mouse, the C-terminal six residues are cleaved, giving rise to a monomeric, 16 kDa heparin-binding bioactive molecule (aa 17-156) (5-7). A shorter form has been described in human (7). The activity seems to be concentrated in the nine aa's preceding the prosegment (aa 148-156). Retention of the prosegment blocks activity (4). The 140 aa mature segment is known to bind to the G-protein coupled receptor termed ChemR23 (5, 7). Binding results in macrophage and immature dendritic cell chemotaxis (5). The distribution of this receptor is limited to immune APCs, and it is assumed that Chemerin is an inflammatory molecule. It is unclear which cells are actually producing Chemerin, but keratinocytes, endothelial cells and osteoclasts are potential candidates (1, 7). Mature mouse Chemerin shares 67%, 84% and 82% aa sequence identity with human, rat and hamster Chemerin, respectively (6). There is apparently cross-species activity for the protein (6).

References:

1. Nagpal, S. *et al.* (1997) *J. Invest. Dermatol.* **109**:91.
2. Storici, P. *et al.* (1996) *Eur. J. Biochem.* **238**:769.
3. Zanetti, M. (2004) *J. Leukoc. Biol.* **75**:39.
4. Wittamer, V. *et al.* (2004) *J. Biol. Chem.* **279**:9956.
5. Wittamer, V. *et al.* (2003) *J. Exp. Med.* **198**:977.
6. Busmann, A. *et al.* (2004) *J. Chromatog. B* **811**:217.
7. Meder, W. *et al.* (2003) *FEBS Lett.* **555**:495.
8. Yokoyama-Kobayashi, M. *et al.* (1999) *Gene* **228**:161.

PRODUCT SPECIFIC NOTICES

This product and/or its use is the subject of European Patent 1 405 083 B1, US Patents 7,332,291; 7,419,658 and 7,842,453 as well as foreign equivalents licensed to R&D Systems Inc. The purchase of this product is intended for research purposes only, not including the screening of compounds for the development of therapeutic and/or diagnostic products. Buyers may require a separate license to the patent rights for applications beyond such research purposes. For information on licensing please contact Euroscreen SA rue Adrienne Bolland n°47 B-6041 Gosselies Belgium. Phone: +32-71-348500, Fax: +32-71-348519, e-mail LicensingBD@euroscreen.com. Attention: Dr. Vincent Lannoy.