

Human IL-18/IL-1F4 Antibody

Monoclonal Mouse IgG₁ Clone # 925008 Catalog Number: MAB2548

DESCRIPTION	
Species Reactivity	Human
Specificity	Detects human IL-18/IL-1F4 in direct ELISAs.
Source	Monoclonal Mouse IgG ₁ Clone # 925008
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	E. coli-derived recombinant human IL-18/IL-1F4 Tyr37-Asp193 Accession # Q14116
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. See Certificate of Analysis for details. *Small pack size (-SP) is supplied either lyophilized or as a 0.2 µm filtered solution in PBS.


APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. General Protocols are available in the Technical Information section on our website.


	Recommended Concentration	Sample
Immunocytochemistry	5-25 μg/mL	See Below
Intracellular Staining by Flow Cytometry	0.25 μg/10 ⁶ cells	See Below
CyTOF-ready	Ready to be labeled with conjugation.	I using established conjugation methods. No BSA or other carrier proteins that could interfere
ELISA	,	ons as an ELISA detection antibody when paired with Rabbit Anti-Human IL-18/IL-1F4 y (Catalog # MAB91243).
	,	nded for assay development on various assay platforms requiring antibody pairs. We recommend 18 DuoSet ELISA Kit (Catalog # DY318-05) for convenient development of a sandwich ELISA oi

DATA

Immunocytochemistry




IL-18/IL-1F4 in Mouse Splenocytes. IL-18/IL-1F4 was detected in immersion fixed mouse splenocytes stimulated with Cal and PMA using Mouse Anti-Human IL-18/IL-1F4 Monoclonal Antibody (Catalog # MAB2548) at 25 µg/mL for 3 hours at room temperature. Cells were stained using the NorthernLights™ 557-conjugated Anti-Mouse IgG Secondary Antibody (red; Catalog # NL007) and counterstained with DAPI (blue). Specific staining was localized to cytoplasm. View our protocol for Fluorescent ICC Staining of Non-adherent Cells



the Human Total IL-18/IL-1F4 Quantikine ELISA Kit (Catalog # DL180) for a complete optimized ELISA.

Detection of IL-18/IL-1F4 in Human THP-1 cell line by Flow Cytometry. THP-1 Human acute monocytic leukemia cell line was stained with Mouse Anti-Human IL-18/IL-1F4 Monoclonal Antibody (Catalog # MAB2548, filled histogram) or isotype control antibody (Catalog # MAB002, open histogram), followed by PE-conjugated Anti-Mouse IgG Secondary Antibody (Catalog # F0102B). To facilitate intracellular staining, cells were fixed with paraformaldehyde and permeabilized with saponin. View our protocol for Staining Intracellular Molecules.



Human IL-18/IL-1F4 ELISA Standard Curve. Recombinant Human IL-18/IL-1F4 protein was serially diluted 2-fold and captured by Rabbit Anti-Human IL-18/IL-1F4 Monoclonal Antibody (Catalog # MAB91243) coated on a Clear Polystyrene Microplate (Catalog # DY990). Mouse Anti-Human IL-18/IL-1F4 Monoclonal Antibody (Catalog # MAB2548) was biotinylated and incubated with the protein captured on the plate. Detection of the standard curve was achieved by incubating Streptavidin-HRP (Catalog # DY998) followed by Substrate Solution (Catalog # DY999) and stopping the enzymatic reaction with Stop Solution (Catalog # DY994).

Rev. 5/10/2018 Page 1 of 2





Human IL-18/IL-1F4 Antibody

Monoclonal Mouse IgG₁ Clone # 925008 Catalog Number: MAB2548

Reconstitution	Reconstitute at 0.5 mg/mL in sterile PBS.	
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below. *Small pack size (-SP) is shipped with polar packs. Upon receipt, store it immediately at -20 to -70 °C	
Stability & Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. 12 months from date of receipt, -20 to -70 °C as supplied. 1 month, 2 to 8 °C under sterile conditions after reconstitution. 6 months, -20 to -70 °C under sterile conditions after reconstitution.	

BACKGROUND

Interleukin-18 (IL-18), also known as IL-1F4 and IFN- γ inducing factor (IGIF), is a member of the IL-1 family of cytokines and is a key molecule in the innate immune response (1). Rhesus IL-18 is synthesized as a 24 kDa proprotein that contains a 36 amino acid (aa) propeptide and a 157 aa mature region (2). Under inflammatory conditions, the propeptide is cleaved by Caspase-1 in the cytoplasm to liberate the mature nonglycosylated 18 kDa monomeric IL-18 (3, 4). Mature rhesus IL-18 shares 96% aa sequence identity with human IL-18 and 60-76% with mouse, rat, canine, feline, and porcine IL-18. IL-18 is secreted by a variety of cell types including macrophages, dendritic cells, and epithelial cells (1, 5). Circulating mature IL-18 is sequestered by soluble IL-18 binding proteins (IL-18 BP) that inhibit IL-18 bioactivity (6). IL-18 interacts with the widely expressed IL-18 Ra which then recruits the signaling subunit IL-18 R β (7, 8). The IL-1 family member IL-177 also binds to IL-18 R β but does not recruit IL-18 R β or induce signaling (9). IL-177 binds IL-18 BP and enhances its neutralizing effect on IL-18 activity (9). IL-18 synergizes with other cytokines to activate NK, Th1, and Th17 cells and to increase the production of IFN- γ (1, 5, 10-12). IL-18 can also promote Th2 cytokine release which reduces the effectiveness of antiviral responses (13, 14). Increased levels of active IL-18 stimulates Th1 and NK cells to target tumor cells, but it can also promote angiogenesis, metastasis, and tumor cell immune evasion (11).

References:

- 1. Arend, W.P. et al. (2008) Immunol. Rev. 223:20.
- 2. Giavedoni, L.D. et al. (2001) J. Interferon Cytokine Res. 21:173.
- 3. Ghayur, T. et al. (1997) Nature 386:619.
- Gu, Y. et al. (1997) Science 275:206.
- 5. Boraschi, D. and C.A. Dinarello (2006) Eur. Cytokine Netw. 17:224.
- Novick, D. et al. (1999) Immunity 10:127.
- 7. Torigoe, K. et al. (1997) J. Biol. Chem. 272:25737.
- 8. Born, T.L. *et al.* (1998) J. Biol. Chem. **273**:29445.
- 9. Bufler, P. et al. (2002) Proc. Natl. Acad. Sci. USA 99:13723.
- 10. Takeda, K. et al. (1998) Immunity 8:383.
- 11. Park, S. et al. (2007) Cell. Mol. Immunol. 4:329.
- 12. Yoshimoto, T. et al. (1998) J. Immunol. 161:3400.
- 13. Hoshino, T. et al. (2001) J. Immunol. 166:7014.
- 14. Iannello, A. et al. (2009) AIDS Rev. 11:115.
- 15. Rabkin, S.W. (2009) Nat. Clin. Pract. Cardiovasc. Med. 6:192.
- 16. Netea, M.G. et al. (2006) Nat. Med. 12:650.

