

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human CD24 in direct ELISAs.
Source	Monoclonal Mouse IgG _{2B} Clone # 528807
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	Chinese hamster ovary cell line CHO-derived recombinant human CD24 Ser27-Gly59 (predicted) Accession # P25063
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. See Certificate of Analysis for details. *Small pack size (-SP) is supplied as a 0.2 µm filtered solution in PBS.


APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

	Recommended Concentration	Sample
Immunohistochemistry	8-25 µg/mL	See Below


DATA

Immunohistochemistry



CD24 in Human Colon. CD24 was detected in immersion fixed paraffin-embedded sections of human colon using 15 µg/mL Human CD24 Monoclonal Antibody (Catalog # MAB5248) overnight at 4 °C. Tissue was stained with the Anti-Mouse HRP-DAB Cell & Tissue Staining Kit (brown; Catalog # CTS002) and counterstained with hematoxylin (blue). View our protocol for [Chromogenic IHC Staining of Paraffin-embedded Tissue Sections](#).

Immunohistochemistry



CD24 in Human Ovarian Cancer Tissue. CD24 was detected in immersion fixed paraffin-embedded sections of human ovarian cancer tissue using Human CD24 Monoclonal Antibody (Catalog # MAB5248) at 15 µg/mL overnight at 4 °C. Tissue was stained using the Anti-Mouse HRP-DAB Cell & Tissue Staining Kit (brown; Catalog # CTS002) and counterstained with hematoxylin (blue). View our protocol for [Chromogenic IHC Staining of Paraffin-embedded Tissue Sections](#).

PREPARATION AND STORAGE

Reconstitution	Reconstitute at 0.5 mg/mL in sterile PBS.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below. *Small pack size (-SP) is shipped with polar packs. Upon receipt, store it immediately at -20 to -70 °C
Stability & Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> ● 12 months from date of receipt, -20 to -70 °C as supplied. ● 1 month, 2 to 8 °C under sterile conditions after reconstitution. ● 6 months, -20 to -70 °C under sterile conditions after reconstitution.

BACKGROUND

CD24, also known as heat stable antigen (HSA), is a heavily and variably glycosylated 35-60 kDa GPI-linked sialoprotein that is expressed on B cells, T cells, keratinocytes, and myofiber synaptic nuclei and is up-regulated in a wide variety of cancers. It is also a marker for exosomes released into the urine and amniotic fluid. CD24 binds to P-Selectin on activated platelets and vascular endothelial cells. CD24 ligation induces tumor cell apoptosis and is important for the deletion of autoreactive T cells. Mature human CD24 shares 30% and 42% aa sequence identity with mouse and rat CD24, respectively.