

# **Human IL-18/IL-1F4 Antibody**

Recombinant Monoclonal Rabbit IgG Clone # 1072H Catalog Number: MAB91244

DESCRIPTION	la contra de la cont
Species Reactivity	Human
Specificity	Detects human IL-18/IL-1F4 in direct ELISAs. In a sandwich immunoassay, with Goat Anti-Human IL-18 BPa Antigen Affinity Purified Polyclonal Antibody (Catalog # AF119) the human IL-18/IL-18 BPa complex is detected.
Source	Recombinant Monoclonal Rabbit IgG Clone # 1072H
Purification	Protein A or G purified from cell culture supernatant
Immunogen	E. coli-derived recombinant human IL-18/IL-1F4 Tyr37-Asp193 Accession # Q14116
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. See Certificate of Analysis for details. *Small pack size (-SP) is supplied either lyophilized or as a 0.2 µm filtered solution in PBS.

#### APPLICATIONS


Please Note: Optimal dilutions should be determined by each laboratory for each application. General Protocols are available in the Technical Information section on our website.

### **ELISA**

This antibody functions as an ELISA capture antibody when paired with Goat Anti-Human IL-18 BPa Antigen Affinity-purified Polyclonal Antibody (Catalog # AF119). This antibody pair detects the human IL-18/IL-18 BPa complex.

This product is intended for assay development on various assay platforms requiring antibody pairs. We recommend the Human IL-18/IL-18 BPa Complex DuoSet ELISA Kit (Catalog # DY8936-05) for convenient development of a sandwich ELISA or the Human IL-18 BPa Quantikine ELISA Kit (Catalog # DBP180) for a complete optimized ELISA.





Human IL-18/IL-18 BPa Complex ELISA Standard Curve. Recombinant Human IL-18/IL-18 BPa complex protein was serially diluted 2-fold and captured by Rabbit Anti-Human IL-18/IL-1F4 Monoclonal Antibody (Catalog # MAB91244) coated on a Clear Polystyrene Microplate (Catalog # DY990). Goat Anti-Human IL-18 BPa Antigen Affinity-purified Polyclonal Antibody (Catalog # AF119) was biotinylated and incubated with the protein captured on the plate. Detection of the standard curve was achieved by incubating Streptavidin-HRP (Catalog # DY998) followed by Substrate Solution (Catalog # DY999) and stopping the enzymatic reaction with Stop Solution (Catalog # DY994).

## PREPARATION AND STORAGE

Reconstitution	Reconstitute at 0.5 mg/mL in steme PBS.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.
	*Small pack size (-SP) is shipped with polar packs. Upon receipt, store it immediately at -20 to -70 °C.

## Stability & Storage

Use a manual defrost freezer and avoid repeated freeze-thaw cycles.

- 12 months from date of receipt, -20 to -70 °C as supplied.
- 1 month, 2 to 8 °C under sterile conditions after reconstitution.
- 6 months, -20 to -70 °C under sterile conditions after reconstitution.

Rev. 1/18/2018 Page 1 of 2





# **Human IL-18/IL-1F4 Antibody**

Recombinant Monoclonal Rabbit IgG Clone # 1072H Catalog Number: MAB91244

### BACKGROUND

Interleukin-18 (IL-18) is a proinflammatory cytokine in the IL-1 family that exerts distinct immune effects depending on the local cytokine environment. It is expressed as a 24 kDa precursor by endothelial and epithelial cells, keratinocytes, γδ T cells, and phagocytes. The precursor is activated intracellularly by Caspase-1 mediated proteolysis to release the 17 kDa mature cytokine. The precursor can also be released by necrotic cells for extracellular cleavage by multiple proteases. IL-18 activation is induced by infection or tissue damage and contributes to disease pathology in chronic inflammation (1-3). IL-18 binds to the widely expressed IL-18 Rα which recruits IL-18 Rβ to form the signaling receptor complex (4, 5). Its bioactivity is negatively regulated by interactions with IL-18 binding proteins and virally encoded IL-18BP homologs (6). In the presence of IL-12 or IL-15, IL-18 enhances anti-viral Th1 immune responses by inducing IFN-γ production and the cytolytic activity of CD8<sup>+</sup> T cells and NK cells (7, 8). In the absence of IL-12 or IL-15, however, IL-18 promotes production of the Th2 cytokines IL-4 and IL-13 by CD4<sup>+</sup> T cells and basophils (9, 10). In the presence of IL-1β or IL-23, IL-18 induces the antigen-independent production of IL-17 by γδ T cells and CD4<sup>+</sup> T cells (11). IL-18 also promotes myeloid dendritic cell maturation and triggers neutrophil respiratory burst (12, 13). In cancer, IL-18 exhibits diverse activities including enhancing anti-tumor immunity, inhibiting or promoting angiogenesis, and promoting tumor cell metastasis (14). Mature human IL-18 shares approximately 63% amino acid sequence identity with mouse and rat IL-18 (15). Alternative splicing in human ovarian cancer generates an isoform that is resistant to Caspase-1 activation (16). A cell surface form can be expressed on M-CSF induced macrophages and released in response to bacterial endotoxin (17).

### References:

- 1. Dinarello, C.A. et al. (2013) Front. Immunol. 4:289.
- 2. Smith, D.E. (2011) J. Leukoc. Biol. 89:383.
- 3. Gu, Y. et al. (1997) Science 275:206.
- 4. Torigoe, K. et al. (1997) J. Biol. Chem. 272:25737.
- 5. Cheung, H. et al. (2005) J. Immunol. 174:5351.
- 6. Novick, D. et al. (1999) Immunity 10:127.
- 7. Fehniger, T.A. et al. (1999) J. Immunol. 162:4511.
- 8. Yoshimoto, T. et al. (1998) J. Immunol. 161:3400.
- 9. Yoshimoto, T. et al. (2000) Nat. Immunol. 1:132.
- Kroeger, K.M. et al. (2009) J. Leukoc. Biol. 86:769.
- 11. Lalor, S.J. et al. (2011) J. Immunol. 186:5738.
- 12. Li, J. et al. (2004) Cell. Immunol. 227:103.
- 13. Elbim, C. et al. (2005) Clin. Diagn. Lab. Immunol. 12:436.
- 14. Fabbi, M. et al. (2015) J. Leukoc. Biol. 97:665.
- 15. Ushio, S. et al. (1996) J. Immunol. 156:4274.
- 16. Gaggero, A. et al. (2004) Oncogene 23:7552.
- 17. Bellora, F. et al. (2012) Eur. J. Immunol. 42:1618.

Rev. 1/18/2018 Page 2 of 2

