



NorthernLights™ Anti-mouse IgG-NL493

ORDERING INFORMATION

Catalog Number: NL009

Lot Number: LZN11

Size: 0.5 mL

Antibody Concentration: 1 mg/mL

Formulation: Phosphate-buffered saline (PBS)
containing < 0.1% v/v sodium azide

Specificity: Mouse IgG (H+L)

Host: Donkey

Storage: 2-8 °C in the dark

Product Description

This donkey anti-mouse IgG (H+L) antibody was purified by mouse IgG (H+L) affinity chromatography and absorbed against human, rabbit, bovine, goat, and Syrian hamster IgG to eliminate cross-reactivity. The purified antibody is then conjugated to fluorochrome NL493*. The spectral characteristics of NL493 are provided in Table 1. For comparison, the spectral characteristics of Fluorescein (FITC) and Alexa Fluor® 488 are also included.

Table 1.

Fluorochrome	Absorption Maximum (nm)	Emission Maximum (nm)
NL493	493	514
Fluorescein, FITC	492	520
Alexa Fluor 488	495	519

Intended Use

For use as a secondary developing reagent in immunofluorescence assays including immunohistochemistry, immunocytochemistry, flow cytometry, and other fluorescent immunoassays with unlabeled primary antibodies of mouse origin.

Reagent Preparation

Centrifuge the antibody solution in a microfuge before use. For most applications, a dilution of 1:200 is recommended. **Optimal dilutions should be determined by each laboratory for each application.**

Storage

Store the reagent at 2-8 °C **in the dark. Do not freeze.** Under these conditions, the product is stable for at least 6 months from the date of receipt.

Precaution

Contains sodium azide as a preservative. Sodium azide may react with lead and copper plumbing to form explosive metal azides. Flush with large volumes of water during disposal. Dispose of azide containing liquids with caution and according to local regulations.

**Additional fluorochrome-labeled donkey anti-mouse IgG antibodies with different spectral characteristics are also available (see Table 2).*

Table 2.

Catalog Number	Product Description	Absorption Maximum (nm)	Emission Maximum (nm)
NL007	Anti-mouse IgG-NL557	557	574
NL008	Anti-mouse IgG-NL637	637	658

For a complete list of NorthernLights products, please visit www.RnDSystems.com.

For research use only. Not for use in human diagnostic, human therapeutic, or human *in vivo* applications.

R&D Systems, Inc.
1-800-343-7475

NorthernLights is a trademark of R&D Systems, Inc.

Alexa Fluor is a registered trademark of Invitrogen, Inc.