

NFκB-p65(Phospho-Ser468) Antibody

Catalog No: #11013



Package Size: #11013-1 50ul #11013-2 100ul #11013-4 25ul


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NFκB-p65(Phospho-Ser468) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB IF
Species Reactivity	Hu Ms Rt
Immunogen Type	Peptide-KLH
Target Name	NFκB-p65
Modification	Phospho-Ser468
Other Names	p65, NFKB3,
Concentration	1.0mg/ml

Images



Immunofluorescence staining of methanol-fixed HeLa cells using NF-κB p65 (phospho-Ser468) antibody (#11013, Red).



Western blot analysis of extracts using NF-κB p65 (phospho-Ser468) antibody (#11013).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.