

Synaptotagmin 1 (Phospho-Ser309) Antibody

Catalog No: #11209

Package Size: #11209-1 50ul #11209-2 100ul #11209-4 25ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	Synaptotagmin 1 (Phospho-Ser309) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	IF
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous level of Synaptotagmin 1 only when phosphorylated at serine 309.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around phosphorylation site of serine 309 (G-L-S(p)-D-P) derived from Human Synaptotagmin 1.
Target Name	Synaptotagmin 1
Modification	Phospho-Ser309
Other Names	SYT1; Sytl; Synaptotagmin I; P65; SYT
Accession No.	Swiss-Prot: P21579NCBI Protein: NP_001129277.1
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 46kd

Immunofluorescence: 1:100~1:200

Images



Immunofluorescence staining of methanol-fixed HeLa cells using Synaptotagmin 1 (Phospho-Ser309) Antibody #11209.

Background

The synaptotagmins are integral membrane proteins of synaptic vesicles thought to serve as Ca(2+) sensors in the process of vesicular trafficking and exocytosis. Calcium binding to synaptotagmin I participates in triggering neurotransmitter release at the synapse

Gustavsson N, et al. Proc Natl Acad Sci U S A. 2008 Mar 11; 105(10):3992-7.

Cnops L, et al. Cereb Cortex. 2008 May; 18(5):1221-31.

Lynch KL, et al. Mol Biol Cell. 2007 Dec; 18(12):4957-68.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.