

## Acetyl-CoA Carboxylase (Phospho-Ser79) Antibody

Catalog No: #11584




Package Size: #11584-1 50ul #11584-2 100ul #11584-4 25ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Acetyl-CoA Carboxylase (Phospho-Ser79) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB
Species Reactivity	Hu
Immunogen Type	Peptide-KLH
Target Name	Acetyl-CoA Carboxylase
Modification	Phospho-Ser79
Other Names	ACC, ACAC, ACC1, ACCA,
Concentration	1.0mg/ml

## Images



Western blot analysis of extracts from 293 cells untreated or treated with AICAR using Acetyl-CoA Carboxylase (Phospho-Ser79) Antibody #11584.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.