

Btk (Phospho-Tyr551) Polyclonal Antibody

Catalog No: #12223

Package Size: #12223-1 50ul #12223-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Btk (Phospho-Tyr551) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Phospho-Btk (Y551) Polyclonal Antibody detects endogenous levels of Btk protein only when phosphorylated at Y551.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Btk around the phosphorylation site of Y551.
Target Name	Btk
Modification	Phospho-Tyr551
Other Names	BTK; AGMX1; ATK; BPK; Tyrosine-protein kinase BTK; Agammaglobulinaemia tyrosine kinase; ATK; B-cell progenitor kinase; BPK; Bruton tyrosine kinase
Accession No.	Swiss-Prot: Q06187NCBI Gene ID: 695
Target Species	human
SDS-PAGE MW	80kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western blotting: 1/500 - 1/2000

ELISA: 1/40000

Not yet tested in other applications

Images



Western Blot analysis of HELA cells using Phospho-Btk (Y551) Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.