

# Ron (Phospho-Ser1394) Polyclonal Antibody

Catalog No: #12239

Package Size: #12239-1 50ul #12239-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Ron (Phospho-Ser1394) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms
Specificity	Phospho-Ron (S1394) Polyclonal Antibody detects endogenous levels of Ron protein only when phosphorylated at S1394.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Ron around the phosphorylation site of S1394.
Target Name	Ron
Modification	Phospho-Ser1394
Other Names	MST1R; PTK8; RON; Macrophage-stimulating protein receptor; MSP receptor; CDw136; Protein-tyrosine kinase 8; p185-Ron; CD antigen CD136
Accession No.	Swiss-Prot: Q04912NCBI Gene ID: 4486
Target Species	human
SDS-PAGE MW	152kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western blotting: 1/500 - 1/2000

ELISA: 1/40000

Not yet tested in other applications

## Images

Western Blot analysis of 3T3 cells using Phospho-Ron (S1394) Polyclonal Antibody



Note: This product is for in vitro research use only and is not intended for use in humans or animals.