

4.1R (Phospho-Tyr660) Polyclonal Antibody

Catalog No: #12262

Package Size: #12262-1 50ul #12262-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	4.1R (Phospho-Tyr660) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms
Specificity	Phospho-4.1R (Y660) Polyclonal Antibody detects endogenous levels of 4.1R protein only when phosphorylated at Y660.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human 4.1R around the phosphorylation site of Y660.
Target Name	4.1R
Modification	Phospho-Tyr660
Other Names	EPB41; E41P; Protein 4.1; P4.1; 4.1R; Band 4.1; EPB4.1
Accession No.	Swiss-Prot: P11171NCBI Gene ID: 2035
Target Species	human
SDS-PAGE MW	60kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western blotting: 1/500 - 1/2000

Immunohistochemistry: 1/100 - 1/300


ELISA: 1/5000

Not yet tested in other applications

Images



Western Blot analysis of JK cells using Phospho-4.1R (Y660) Polyclonal Antibody



Western Blot analysis of Jurkat cells using Phospho-4.1R (Y660) Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.