

Flg (Phospho-Tyr654) Polyclonal Antibody

Catalog No: #12267



Package Size: #12267-1 50ul #12267-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Flg (Phospho-Tyr654) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	Phospho-Flg (Y654) Polyclonal Antibody detects endogenous levels of Flg protein only when phosphorylated at Y654.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Flg around the phosphorylation site of Y654.
Target Name	Flg
Modification	Phospho-Tyr654
Other Names	FGFR1; BFGFR; CEK; FGFBR; FLG; FLT2; HBGFR; Fibroblast growth factor receptor 1; FGFR-1; Basic fibroblast growth factor receptor 1; BFGFR; bFGF-R-1; Fms-like tyrosine kinase 2; FLT-2; N-sam; Proto-oncogene c-Fgr; CD antigen CD331
Accession No.	Swiss-Prot: P11362NCBI Gene ID: 2260
Target Species	human
SDS-PAGE MW	120kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1/500 - 1/2000


Immunofluorescence: 1/200 - 1/1000

ELISA: 1/20000

Not yet tested in other applications

Images

Western Blot analysis of HeLa cells using Phospho-Flg (Y654)
Polyclonal Antibody



Note: This product is for in vitro research use only and is not intended for use in humans or animals.