

CD5 (Phospho-Tyr453) Polyclonal Antibody

Catalog No: #12269



Package Size: #12269-1 50ul #12269-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD5 (Phospho-Tyr453) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	Phospho-CD5 (Y453) Polyclonal Antibody detects endogenous levels of CD5 protein only when phosphorylated at Y453.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human CD5 around the phosphorylation site of Y453.
Target Name	CD5
Modification	Phospho-Tyr453
Other Names	CD5; LEU1; T-cell surface glycoprotein CD5; Lymphocyte antigen T1/Leu-1; CD antigen CD5
Accession No.	Swiss-Prot: P06127NCBI Gene ID: 921
Target Species	human
SDS-PAGE MW	40kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1/500 - 1/2000

Immunohistochemistry: 1/100 - 1/300


Immunofluorescence: 1/200 - 1/1000

ELISA: 1/20000

Not yet tested in other applications

Images

Western Blot analysis of 293 cells using Phospho-CD5 (Y453)
Polyclonal Antibody



Note: This product is for in vitro research use only and is not intended for use in humans or animals.