

Rb (Phospho-Thr826) Polyclonal Antibody

Catalog No: #12274



Package Size: #12274-1 50ul #12274-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rb (Phospho-Thr826) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	Phospho-Rb (T826) Polyclonal Antibody detects endogenous levels of Rb protein only when phosphorylated at T826.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Rb around the phosphorylation site of T826.
Target Name	Rb
Modification	Phospho-Thr826
Other Names	RB1; Retinoblastoma-associated protein; p105-Rb; pRb; Rb; pp110
Accession No.	Swiss-Prot: P06400NCBI Gene ID: 5925
Target Species	human
SDS-PAGE MW	110kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1/500 - 1/2000


Immunohistochemistry: 1/100 - 1/300

Immunofluorescence: 1/200 - 1/1000

ELISA: 1/10000

Not yet tested in other applications

Images



Western blot analysis of RAT-muscle using p-Rb (T826) antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.