

## FADD (Phospho-Ser194) Polyclonal Antibody

Catalog No: #12299



Package Size: #12299-1 50ul #12299-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	FADD (Phospho-Ser194) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu
Specificity	Phospho-FADD (S194) Polyclonal Antibody detects endogenous levels of FADD protein only when phosphorylated at S194.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human FADD around the phosphorylation site of S194.
Target Name	FADD
Modification	Phospho-Ser194
Other Names	FADD; MORT1; GIG3; Protein FADD; FAS-associated death domain protein; FAS-associating death domain-containing protein; Growth-inhibiting gene 3 protein; Mediator of receptor induced toxicity
Accession No.	Swiss-Prot: Q13158NCBI Gene ID: 8772
Target Species	human
SDS-PAGE MW	28kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western blotting: 1/500 - 1/2000

Immunohistochemistry: 1/100 - 1/300

ELISA: 1/5000

Not yet tested in other applications

## Images



Western blot analysis of MOUSE-KIDNEY MOUSE-BRAIN  
MOUSE-SPLEEN MOUSE-LUNG using p-FADD (S194)  
antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.