

Vimentin (Phospho-Ser83) Polyclonal Antibody

Catalog No: #12333



Package Size: #12333-1 50ul #12333-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vimentin (Phospho-Ser83) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	Phospho-Vimentin (S83) Polyclonal Antibody detects endogenous levels of Vimentin protein only when phosphorylated at S83.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Vimentin around the phosphorylation site of S83.
Target Name	Vimentin
Modification	Phospho-Ser83
Other Names	VIM; Vimentin
Accession No.	Swiss-Prot: P08670NCBI Gene ID: 7431
Target Species	human
SDS-PAGE MW	54kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western blotting: 1/500 - 1/2000

Immunohistochemistry: 1/100 - 1/300

ELISA: 1/10000

Not yet tested in other applications

Images



Western Blot analysis of 293 cells using Phospho-Vimentin (S83) Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.