

β²-Actin Antibody

Catalog No: #21338



Package Size: #21338-1 50ul #21338-2 100ul #21338-4 25ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	β ² -Actin Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB
Species Reactivity	Hu Ms Rt
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around aa.6-10(A-A-L-V-V)derived from Human beta-actin.
Target Name	β ² -Actin
Other Names	BRWS1, PS1TP5BP1,
Concentration	1.0mg/ml

Application Details

Predict MW: 43kd

Western Blot: 1:2000-1:3000

Images



Western blot analysis of extracts from 293, HeLa and 3T3 cells using β-actin Antibody #21338.

Background

Actins are highly conserved proteins that are involved in various types of cell motility and are ubiquitously expressed in all eukaryotic cells.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.