

## RNA silencing suppressor Antibody

Catalog No: #21411



Package Size: #21411-1 50ul #21411-2 100ul #21411-4 25ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	RNA silencing suppressor Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	WB
Species Reactivity	Virus
Specificity	The antibody detects endogenous level of total RNA silencing suppressor protein.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around aa.9~13 (D-A-R-E-Q) derived from Cucumber necrosis virus RNA silencing suppressor.
Target Name	RNA silencing suppressor
Accession No.	Swiss-Prot: P15184 NCBI Protein: NP_040957.1
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

## Application Details

Predicted MW: 20kd

Western blotting: 1:500~1:1000

## Background

plays a role in initiation of gene silencing at the early stage of viral infection.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.