

syk(Ab-323) Antibody

Catalog No: #21546



Package Size: #21546-1 50ul #21546-2 100ul #21546-4 25ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	syk(Ab-323) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	WB IHC IF
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total syk protein.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around aa. 321~325 (N-P-Y-E-P) derived from Human syk.
Target Name	syk
Other Names	Spleen tyrosine kinase
Accession No.	Swiss-Prot: P43405NCBI Protein: NP_001128524.1
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 72kd
Western blotting: 1:500~1:1000
Immunohistochemistry: 1:50~1:100
Immunofluorescence: 1:100~1:200


Images



Western blot analysis of extracts from JK cells using syk(Ab-323) Antibody #21546.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using syk(Ab-323) Antibody #21546(left) or the same antibody preincubated with blocking peptide(right).



Immunofluorescence staining of methanol-fixed HeLa cells using syk(Ab-323) Antibody #21546.

Background

Positive effector of BCR-stimulated responses. Couples the B-cell antigen receptor (BCR) to the mobilization of calcium ion either through a phosphoinositide 3-kinase-dependent pathway, when not phosphorylated on tyrosines of the linker region, or through a phospholipase C-gamma-dependent pathway, when phosphorylated on Tyr-348 and Tyr-352. Thus the differential phosphorylation of Syk can determine the pathway by which BCR is coupled to the regulation of intracellular calcium ion

Zhang, J. et al. (2000) J. Biol. Chem. 275, 35442-35447.

Turner, M. et al. (2000) Immunol. Today 21, 148-154.

Decker, M. et al. (1998) J. Biol. Chem. 273, 8867-8874.

Law, C.L. et al. (1996) Mol. Cell. Biol. 16, 1305-1315.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.