

Puma Antibody

Catalog No: #21688



Package Size: #21688-1 50ul #21688-2 100ul #21688-4 25ul


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Puma Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB
Species Reactivity	Hu
Immunogen Type	Peptide-KLH
Target Name	Puma
Other Names	BBC3, JFY1, JFY-1,
Concentration	1.0mg/ml

Images



Western blot analysis of extracts from 293 and MCF-7 cells using Puma Antibody #21688.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.