

SGTA antibody

Catalog No: #22393

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


| | |
|-----------------------|--|
| Product Name | SGTA antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | Purified by antigen-affinity chromatography. |
| Applications | WB |
| Species Reactivity | Hu |
| Immunogen Type | Recombinant protein |
| Immunogen Description | Recombinant protein fragment contain a sequence corresponding to a region within amino acids 6 and 313 of Human SGTA |
| Target Name | SGTA |
| Accession No. | NCBI Gene ID: 6449NCBI mRNA#: NM_003021NCBI Protein#: NP_003012 |
| Formulation | Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative. |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |

Application Details

Predicted MW: 34kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)

A: Hep G2

B: Molt-4

10% SDS PAGE

Primary antibody diluted at 1: 5000

Background

This gene encodes a protein which is capable of interacting with the major nonstructural protein of parvovirus H-1 and 70-kDa heat shock cognate protein; however, its function is not known. Since this transcript is expressed ubiquitously in various tissues, this protein may serve a housekeeping function. [provided by RefSeq]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.