

TFEB antibody

Catalog No: #22710

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	TFEB antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 1 and 14 of Human TFEB
Target Name	TFEB
Other Names	AlphaTFEB; TCFEB
Accession No.	NCBI Gene ID: 7942NCBI mRNA#: NM_007162NCBI Protein#: NP_009093
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 53kd

Western blotting: 1:500-1:3000

Images



Sample(30 ug whole cell lysate)
A: 293T
10% SDS PAGE
Primary antibody diluted at 1: 1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.