

SARS Spike Antibody

Catalog No: #24216

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SARS Spike Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	E
Species Reactivity	Virus
Immunogen Type	Peptide
Immunogen Description	Raised against a synthetic peptide corresponding to amino acids at the amino-terminus of. The SARS Spike glycoprotein.
Target Name	SARS Spike
Other Names	SARS
Accession No.	P59594
Formulation	Supplied in PBS containing 0.02% sodium azide.
Storage	Can be stored at -20°C, stable for one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Background

A novel coronavirus has recently been identified as the causative agent of SARS (Severe Acute Respiratory Syndrome). Coronaviruses are a major cause of upper respiratory diseases in humans. The genomes of these viruses are positive-stranded RNA approximately 27-31kb in length. SARS infection can be mediated by the binding of the viral spike protein, a glycosylated 139 kDa protein and the major surface antigen of the virus, to the angiotensin-converting enzyme 2 (ACE2) on target cells. This binding can be blocked by a soluble form of ACE2.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.