

Rift Valley Fever Virus Antibody

Catalog No: #24652

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rift Valley Fever Virus Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	E
Species Reactivity	Virus
Immunogen Type	Peptide
Immunogen Description	Raised against a 20 amino acid peptide near the center of the RVF virus.
Target Name	Rift Valley Fever Virus
Other Names	RVF virus, RVFV
Accession No.	P03518
Formulation	Supplied in PBS containing 0.02% sodium azide.
Storage	Can be stored at -20°C, stable for one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Background

Rift Valley Fever (RVF) virus is an arthropod-borne virus endemic to Africa that infects humans and animals that is transmitted predominantly by mosquitoes. During human infections, symptoms can range from benign fever to severe encephalitis and fatal hepatitis with hemorrhagic fever. The Bunyaviridae family of viruses to which the RVF virus belongs are spherical enveloped viruses with a tripartite RNA genome of negative or ambisense polarity. The three segments are referred to as the L, M, and S segments. The L and M segments are negative polarity and code for the L-dependent RNA polymerase and glycoprotein precursor respectively. The S segment is of ambisense polarity and encodes the nucleoprotein and non-structural proteins. This RVF virus antibody was derived from a peptide sequence near the center of the polyprotein precursor translated from the M segment. It will therefore detect both the precursor and the Glycoprotein G1.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.