

NUMB Antibody


Catalog No: #25307

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NUMB Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	E WB
Species Reactivity	Hu Ms Rt
Specificity	At least four isoforms of NUMB are known to exist; this antibody will detect all of them.
Immunogen Type	Peptide
Immunogen Description	Raised against a 17 amino acid peptide near the carboxy terminus of human NUMB.
Target Name	NUMB
Other Names	m-NUMB, S171
Accession No.	NP_1005743
Formulation	Supplied in PBS containing 0.02% sodium azide.
Storage	Can be stored at -20°C, stable for one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Images



Western blot analysis of NUMB in mouse lung tissue lysate with NUMB antibody at (A) 0.25 and (B) 0.5 ug/mL.

Background

NUMB, the mammalian homolog to the *Drosophila* asymmetric cell fate determinant NUMB, is thought to share several features molecular mechanisms in mammalian cells, generating asymmetric cell divisions during neurogenesis in vertebrate development as well as in hematopoietic stem cells. NUMB has been shown to inhibit Notch signaling, and is itself regulated by ubiquitylation by MDM2. NUMB has also been shown to help activate the tumor suppressor p53, suggesting that loss of NUMB in cancerous cells would not only activate the potential oncogene Notch, but diminish the tumor suppressing effect of p53.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.