

BRUCE Antibody

Catalog No: #25439

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	BRUCE Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	E IHC
Species Reactivity	Hu
Specificity	BRUCE antibody is human specific. Multiple isoforms of BRUCE are known to exist.
Immunogen Type	Peptide
Immunogen Description	17 amino acid peptide near the carboxy terminus of human BRUCE.
Target Name	BRUCE
Other Names	Baculoviral IAP repeat containing 6, BIRC6, APOLLON
Accession No.	NP_057336
Formulation	PBS containing 0.02% sodium azide.
Storage	Store at 4°C for three months and -20°C, stable for up to one year.

Application Details

Predicted MW: 534kd

Images



Immunohistochemistry of BRUCE in human small intestine tissue with BRUCE antibody at 5 ug/ml.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.