

Pst1 Antibody

Catalog No: #25443

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Pst1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	E WB
Species Reactivity	Yeast
Immunogen Type	Peptide
Immunogen Description	synthetic peptide from the amino terminus of yeast Pst1.
Target Name	Pst1
Other Names	Protoplast secreted protein 1
Accession No.	Q12355
Formulation	PBS containing 0.02% sodium azide.
Storage	Store at 4°C for three months and -20°C, stable for up to one year.

Application Details

Predicted MW: 49kd

Images



Western blot analysis of TM yeast Pst1 protein (50 ng) with Pst1 antibody at 1 ug/ml.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.