

NCLN Antibody


Catalog No: #25471

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NCLN Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	E WB
Species Reactivity	Hu Ms Rt
Immunogen Type	Peptide
Immunogen Description	NCLN antibody was raised against an 18 amino acid peptide near the carboxy terminus of human NCLN.
Target Name	NCLN
Other Names	Nicalin, NET59
Accession No.	NP_064555
Formulation	PBS containing 0.02% sodium azide.
Storage	Store at 4°C for three months and -20°C, stable for up to one year.

Images



Western blot analysis of NCLN in mouse heart tissue lysate with NCLN antibody at 0.5 ug/mL in (A) the absence and (B) the presence of blocking peptide.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.