

NIBRIN Antibody


Catalog No: #25487

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NIBRIN Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	E WB
Species Reactivity	Hu Ms
Specificity	Two alternatively spliced transcript isoforms of NIBRIN are known to exist.
Immunogen Type	Peptide
Immunogen Description	15 amino acid peptide near the carboxy terminus of human NIBRIN.
Target Name	NIBRIN
Other Names	Nijmegen breakage syndrome protein 1, NBN, AT-V1, AT-V2, ATV, NBS, NBS1, P95, Cell cycle regulatory protein p95
Accession No.	NP_002476
Formulation	PBS containing 0.02% sodium azide.
Storage	Store at 4°C for three months and -20°C, stable for up to one year.

Images



Western blot analysis of NIBRIN in mouse lung tissue lysate with NIBRIN antibody at 1 ug/ml.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.