

PTPRD Antibody


Catalog No: #25495

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PTPRD Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	E WB
Species Reactivity	Hu
Specificity	PTPRD antibody is human specific. At least three alternatively spliced transcript variants encoding distinct isoforms have been observed.
Immunogen Type	Peptide
Immunogen Description	PTPRD antibody was raised against an 18 amino acid peptide near the amino terminus of human PTPRD.
Target Name	PTPRD
Other Names	Protein tyrosine phosphatase receptor type D, PTPD, RPTPDELTA, R-PTP-delta
Accession No.	NP_569075
Formulation	PBS containing 0.02% sodium azide.
Storage	Store at 4°C for three months and -20°C, stable for up to one year.

Images



Western blot analysis of PTPRD in HeLa cell lysate with PTPRD antibody at (A) 1 and (B) 2 ug/ml.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.