

SELK Antibody

Catalog No: #25519

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SELK Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	E WB
Species Reactivity	Hu Ms Rt
Immunogen Type	Peptide
Immunogen Description	SELK antibody was raised against an 18 amino acid peptide near the carboxy terminus of human SELK.
Target Name	SELK
Other Names	Selenoprotein K, HSPC030, HSPC297
Accession No.	NP_067060
Formulation	PBS containing 0.02% sodium azide.
Storage	Store at 4°C for three months and -20°C, stable for up to one year.

Application Details

Predicted MW: 10kd

Images



Western blot analysis of SELK in A20 cell lysate with SELK antibody at (A) 1 and (B) 2 ug/ml.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.