

ANG Antibody

Catalog No: #25522

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ANG Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	E WB
Species Reactivity	Hu Ms Rt
Specificity	Two alternatively spliced transcript variants have been observed.
Immunogen Type	Peptide
Immunogen Description	15 amino acid peptide near the center of human ANG.
Target Name	ANG
Other Names	Angiogenin ribonuclease, RNase A family 5, ALS9, angiogenin ribonuclease 5, HEL168, RNASE4, RNASE5
Accession No.	NP_001136
Formulation	PBS containing 0.02% sodium azide.
Storage	Store at 4°C for three months and -20°C, stable for up to one year.

Application Details

Predicted MW: 16kd

Images



Western blot analysis of ANG in rat liver tissue lysate with Ang antibody at (A) 1 and (B) 2 ug/ml.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.