

## EED Antibody

Catalog No: #25531

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	EED Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	E WB
Species Reactivity	Hu Ms Rt
Specificity	Four alternatively spliced transcript variants have been observed.
Immunogen Type	Peptide
Immunogen Description	EED antibody was raised against an 18 amino acid peptide near the amino terminus of human EED.
Target Name	EED
Other Names	Embryonic ectoderm development, HEED, Polycomb protein EED, WAIT1, WAIT-1, WD protein associating with integrin cytoplasmic tails 1
Accession No.	NP_003788
Formulation	PBS containing 0.02% sodium azide.
Storage	Store at 4°C for three months and -20°C, stable for up to one year.

## Application Details

Predicted MW: 49kd

## Images



Western blot analysis of EED in human heart tissue lysate with Eed antibody at 1 ug/ml.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.