

JMJD6(N-term) Monoclonal Antibody

Catalog No: #27005

Package Size: #27005-1 50ul #27005-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	JMJD6(N-term) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Purified by antigen-affinity chromatography.
Applications	WB ICC
Species Reactivity	Hu D Ms Rt
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids 11-248aa of human JMJD6(N-terminus).
Target Name	JMJD6
Other Names	bifunctional arginine demethylase and lysyl-hydroxylase JMJD6 <input type="checkbox"/> lysyl-hydroxylase JMJD6 <input type="checkbox"/> phosphatidylserine receptor <input type="checkbox"/> jmjC domain-containing protein 6 <input type="checkbox"/> histone arginine demethylase JMJD6 <input type="checkbox"/> peptide-lysine 5-dioxygenase JMJD6 <input type="checkbox"/> jumonji doma
Accession No.	UniProt:Q6NYC1 NCBI Gene ID:23210
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

Application Details

Predicted MW: 62kd


Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected JMJD6 on 40ug of cell lysates

1:200 dilution from a previous lot detected JMJD6 in Hela



Note: This product is for in vitro research use only and is not intended for use in humans or animals.