

# Ikaros(C-term) Monoclonal Antibody

Catalog No: #27062



Package Size: #27062-1 50ul #27062-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	Ikaros(C-term) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu Ms D
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids 249-415aa of human IKAROS(C-terminus).
Target Name	Ikaros
Other Names	DNA-binding protein Ikaros □ CLL-associated antigen KW-6 □ Ikaros (zinc finger protein) □ ikaros family zinc finger protein 1 □ lymphoid transcription factor LyF-1 □ zinc finger protein, subfamily 1A, 1 (Ikaros
Accession No.	UniProt:Q13422 NCBI Gene ID:10320
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

## Application Details

Predicted MW: 58kd

Western blotting: 1:1000

## Images



1:1000 dilution of this antibody detected IKAROS on 40ug of cell lysates

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.