

Integrin beta 3(N-terminus) Monoclonal Antibody

Catalog No: #27094



Package Size: #27094-1 50ul #27094-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	Integrin beta 3(N-terminus) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	3C12-G2-E4-H9
Isotype	IgG1
Purification	Affinity purified
Applications	E WB
Species Reactivity	Hu D Rb Rt Ms
Immunogen Type	Recombinant protein
Immunogen Description	Purified recombinant human ITGB3 protein fragments expressed in E.coli.
Target Name	Integrin beta3
Other Names	GT, CD61, GP3A, BDPLT2, GPIIIa
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol

Application Details

Predicted MW: 120kd

Western blotting: 1:1000

Images



Western blot detection of ITGB3 in Leukocyte whole cell lysate(1:1000 diluted).

Predicted band size:120KDa.

Observed band size:120KDa.

Background

Integrin alpha-V/beta-3 is a receptor for cytotactin, fibronectin, laminin, matrix metalloproteinase-2, osteopontin, osteomodulin, prothrombin, thrombospondin, vitronectin and von Willebrand factor. Integrin alpha-IIb/beta-3 is a receptor for fibronectin.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.