

# NUP98 Monoclonal Antibody

Catalog No: #27111



Package Size: #27111-1 50ul #27111-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	NUP98 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	3B8-D7-H10
Isotype	IgG1
Purification	Affinity purified
Applications	E WB
Species Reactivity	Hu D Ms
Immunogen Type	Recombinant protein
Immunogen Description	Purified recombinant human NUP98 protein fragments expressed in E.coli.
Target Name	NUP98
Other Names	ADIR2, NUP96, NUP196
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol

## Application Details

Predicted MW: 98kd


Western blotting: 1:1000

## Images



Western blot detection of NUP98 in HeLa

Immunoprecipitation analysis of HeLa cell lysates using NUP98 antibody



## Background

Nup98 and Nup96 play a role in the bidirectional transport across the nucleoporin complex (NPC). The repeat domain in Nup98 has a direct role in the transport.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.